




HYDRO MASTER

10 IN 1 Hydra Dermabrasion

- Cleanse + Extract + Hydrate
- Boost Your Skin. Boost Your Confidence.
- All-In-One Analyzing Platform for All Kind of Skin.





FEATURE

1. The treatment leaves skin clear, radiant, and hydrated.
2. 4 big bottles for skin solutions, save your time to change bottles frequently.
3. Noninvasive no pain, no injury.
4. Full Range of Liquid: It can use any liquid (such as makeup water, energy water).
5. A complete effective, sterile method of treatment anti inflammatory.

FUNCTION



Shrink Pores



Deep Cleaning



Skin Firming



Blackheads Removing



Skin Care



Dark Circles Removal



Anti-Aging



Brighten Skin Tone






10 Handles



BUBBLE CLEANING

The three major regions and the division of labor are clear, and the effect can be tripled at one time.

Vacuum, with negative pressure sucking and taking away all dirty from the face



Soft and hard replacement heads, targeted cleansing for different skin types



This new generation of skin cleansing device can produce nano-level dense foam that penetrates deep into the pores, solving the problem of clogged pores from the root, leaving the skin clean, fair and soft.

The Oxygen Bubble device is a handheld, ergonomic tool with a grey and white body. It features a power button on the side and a white, dome-shaped nozzle at the tip. The nozzle is shown detached and lying on the surface below the device, revealing a circular opening with several small, raised bumps. The device is set against a dark blue background with glowing blue light trails.

OXYGEN BUBBLE

Activates cells for enhanced absorption, injecting the product deep into the skin; achieving whitening and brightening effects on dull skin.

The Skin Scrubber device is a handheld, ergonomic tool with a grey and white body. It features a power button on the side and a long, thin, curved metal tip at the end. The device is set against a dark blue background with glowing blue light trails.

SKIN SCRUBBER

The skin scrubber uses high-frequency vibrations to soften skin cells and deeply cleanse impurities.
(Improve acne scars, minimize pores, and remove blackheads)

HOT & COLD HAMMER



Cooling function:

Tightens pores, calms and repairs damaged skin, and reduces redness

Thermal function:

Increases skin surface temperature, accelerates blood circulation, and increases product absorption rate.

BIPOLAR HR



High-frequency vibrations create heat through friction within the skin's molecules, thereby tightening the skin and removing wrinkles.



ULTRASONIC MESOTHERAPY

Ultrasonic wave + radio frequency automatic liquid introduction technology allows the product absorption rate to reach 90%, refining pores and improving skin condition



DUAL-CIRCLE THERMAL WAVE

Using 0.5-1M variable frequency heat energy, it can act on different depths of the dermis to achieve the effect of reducing wrinkles and tightening the skin.

DIAMOND DERMABRASION



Removes aging cuticles from the outer layer of the skin to promote skin renewal. The strong suction can effectively absorb oil, dirt, acne, and blackheads.

PLASMA PEN



Use Plasma function can quickly remove and sweep away spots without bleeding, improving skin radiance



**OXYGEN
INJECTION**

High-pressure oxygen injection uses small particles to penetrate deep into the skin, enhancing moisture retention and allowing for easier absorption by shrinking pores.



12.1 INCH SCREEN

- Intelligent operating system
- Easy to operate

BEFORE & AFTER



MACHINE PARAMETERS



VOLTAGE:	100V-120V 50Hz / 220V-240V 60Hz
POWER:	180W ≥
HR FREQUENCY:	0.5-1 Mhz
COLD AND HOT SCOPE :	Cold 0-6°C / Hot 36°C- 46°C
ULTRASOUND:	1.1Mhz
PROBE:	6 soft heads, 18 hard heads
CLEAN WATER MILL:	Positive pressure 2-3Kg, Negative pressure 96Kpa, Flow ≤ 30L
DISPLAY:	12.1 inches
SOLUTION BOTTLE:	A.Deep layer clearing / B.Whitening and moisturizing / C. Deep cleaning / D. Instrument piping cleaning